

Titel **Column screen 123 release**

K5010



Analysis	Requirement	Result
Antigen control	Various antigens have been tested and found conform the enclosed antigram of the panel Tested in: Ortho BioVue System (anti-IgG, anti-C3d; (polyspecific) and in BIO-RAD ID-Card (LISS/Coombs)	Conform requirements
Hematocrit	0.8 ± 0.2%	Conform requirements

All the donor units used in the production of this reagent have been tested for the markers of specific disease transmitting agents in accordance with current EU guidelines to GMP and found to be non reactive.

Date

: 21 March 2024

Approved by QA/QC



VERTROUWELIJK

Use copy according to REA-SOP-0002.

Printed on: 12 Mar 2024 at: 02:07:07 pm ; Printed by M.DEMOL