

# Product specific information



Sanquin Reagents B.V.  
Plesmanlaan 125  
1066 CX Amsterdam  
The Netherlands

Phone: +31 20 5123599  
Fax: +31 20 5123570  
reagents@sanquin.nl  
www.sanquin.org/reagents

**PeliCluster Anti-Granzyme B**

**REF M9174.1**

**REF M9174.5**

343\_v01 08/2024 (en)

*For research use only*

## Application

The monoclonal Antibody (Clone GB11)<sup>1-3</sup> is derived from a single cell clone, ensuring high specificity and uniformity in recognizing a single epitope on the Granzyme B protein. Granzyme B is a 32 kDa serine protease. Granzyme B is abundantly stored in the granules of cytotoxic T lymphocytes and NK cells. Low levels of expression have been reported in granulocytes, B cells, and activated dendritic cells. Granzyme B is crucial for rapid induction of cell death and apoptosis through interaction with mannose-6-phosphate receptor.

## Characteristics

Type: Monoclonal Antibody

Isotype: Mouse IgG1, K

Clone: CLB-GB11

## Composition

Product: PeliCluster Anti-Granzyme B

Buffer: 20 mM Tris, 150 mM NaCl, pH 8.0

Concentration: 2 mg/mL

Preservative: None

## Precautions

For Research Use Only. Not for use in diagnostic procedures. Keep only in original packaging. Product should be stored upright. Leaking or damaged vials must not be used. The reagent cannot be assumed to be free from infectious agents. Care must be taken in the use and disposal of each container and its contents. Waste-disposal, after completion of the test, should be performed according to your laboratory regulations. CoA is available through the website.

## Specimen Collection and Handling

Source: Tissue culture

Purification method: Protein A chromatography

Specificity: 8-25% positive with lymphocytes (cytoplasmatic)

Storage: See Label




Figure 1. This figure presents the flow cytometry analysis of Granzyme B, conducted with FITC. The analysis was performed on a BD Biosciences FACS Canto II system to assess the expression levels of Granzyme B.

**Disclaimer**

FOR RESEARCH USE ONLY, See attachment. Sanquin Reagents B.V. declines all responsibility arising from any deviation thereof

1. Kummer, J. A. et al. Production and Characterization of Monoclonal Antibodies Raised against Recombinant Human Granzymes A and B and Showing Cross Reactions with the Natural Proteins \*. Journal of immunological Methods vol. 163 (1993).
2. Wever, P. C. et al. The CD8 + Granzyme B + T-Cell Subset in Peripheral Blood from Healthy Individuals Contains Activated and Apoptosis-Prone Cells. Immunology vol. 93 (1998).
3. Spaeny-Dekking, E. H. A. et al. Extracellular Granzymes A and B in Humans: Detection of Native Species During CTL Responses In Vitro and In Vivo. The Journal of Immunology 160, 3610–3616 (1998).